

4K MEDICAL GRADE MONITORS

55", 32" & 27" 4K Monitor

3840 X 21609 INPUT +

BT.2020 COMPATIBLE FOR BETTER COLOUR ACCURACY NEARER TO HUMAN EYES +

MULTIMODALITY QUADVIEW POP THROUGH 4DIFFERENT INPUTS +



ADDITIONAL FEATURES

- + Available in 27", 32" & 55" Monitors
- + PIP modes
- + Input - 4K HDMI, 3G SDI, DVI-D, Display port
- + Advanced image processing algorithms
- + Support high dynamic range
- + Picture in picture & Picture by picture

SIZE	27"	32"	55"
Display Technology	TFT/IPS, LED Backlight	TFT/IPS, LED Backlight	TFT/IPS, LED Backlight
Colour Gammut	BT 2020	BT 2020	BT 2020
Contrast ratio	1500:1	1300:1	1100:1
Resolution 4K	3840 X 2160	3840 X 2160	3840 X 2160
4K Inputs	1 X DP 1.2, 1 X HDMI 2.0	1 X DP 1.2, 1 X HDMI 2.0	1 X DP 1.1, 2 X HDMI 2.0
Full HD inputs	1 X DP, 1X3G SDI, 1 X DVI	1 X DP, 1X3G SDI, 1 X DVI	1 X 3G SDI, 1 X DVI
Output format	1 X DVI, 1 X 3G SDI	1 X DVI, 1 X 3G SDI	1 X DVI, 1 X 3G SDI
PIP / PBP	PIP / PBP	PIP / PBP	PIP / PBP
Image rotation & flip	Available	Available	Available
Dimensions (WxHxD)	654 X 411 X 58 mm	883 X 599 X 200 mm	1259 X 733 X 87 mm
Certification	CE, UL-US, BIS	CE, UL-US, BIS	CE, UL-US, BIS
Viewing Angle	178°(V) / 178°(H)	178°(V) / 178°(H)	178°(V) / 178°(H)
Brightness	800 cd/m ²	800 cd/m ²	500 cd/m ²
Weight	7.7 Kg	13.5 Kg	24 Kg
Art no :	M02D. 274K	M02D. 324K	M02D. 554K