

Sinuscopes HD

FEATURES

- ⌘ Multicoated optical system
- ⌘ High Magnification HD optics
- ⌘ Wide Angle optics for wider view in restricted space
- ⌘ Optimum colour rendering
- ⌘ High contrast image even at periphery
- ⌘ Laser welded joints for longer life
- ⌘ Sapphire tips for scratch resistance
- ⌘ Glass cone to optimise light input
- ⌘ Nitrogen filled to avoid internal fogging

4.0 mm



4 mm Sinuscope	Art.No.
----------------	---------

0°	2040.1800H
----	------------

30°	2040.1830H
-----	------------

45°	2040.1845H
-----	------------

70°	2040.1870H
-----	------------

WL	185mm
----	-------

3.0 mm



3 mm Sinuscope	Art.No.
----------------	---------

0°	2030.1500S
----	------------

30°	2030.1530S
-----	------------

WL	150 mm
----	--------